

Assembling Your Oto

Please follow the instructions below to help you assemble and use your Oto for Clinicians.

Set Up Instructions



1

Download Mobile App

CellScope Companion App, can be downloaded on the Apple Store. Once downloaded create an account.

2

Attach Compatible Case

Depending on the device you are using (iPhone or iPod Touch) you will need to attach it's compatible case, To attach, slide the case onto your device. Cases include; iPhone 5/5s full case, iPhone 6/6s/7 half case. *iPod half case sold separately.

3

Connect Oto Attachment

Slide the small white Oto attachment into the camera opening on the case and slide into place. Make sure it is securely in place (you should hear a click.) Oto attaches likewise with all compatible cases listed.

4

Attach Specula Tip

When you're ready to perform an ear exam, place a tip on the end of the Oto attachment and twist lightly to firmly attach it.

5

Launch the App

Launch the mobile application and login using your username and password.

6

Begin the Exam

Begin an exam and follow the instructions within the app to store or share the image or video.

Mobile App Operating Instructions

Start the Ear Exam

To start, select Left or Right to indicate the ear that you will be recording first.

Start the Recording

Insert Tip

Gently insert the appropriately sized tip into the patient's ear canal and locate the Tympanic Membrane.

Capture View

Once you have captured a view of the Tympanic Membrane on the device screen, end the recording. Note: there is no required video length.

Repeat for Other Ear

Repeat steps 1-4 for the other ear.

Additional Information

Follow the instructions within the app to add additional information and save the exam.

Detach the Oto

When the exam is complete, detach the Oto from the case and store safely, until your next ear exam.

1



2

3



4



5

6



7