

Carbon Fiber Butterfly Boards

Versatile positioning devices for breast, lung and thorax treatments.

The Bionix Butterfly Boards are lightweight and durable devices that comfortably support a patient's arms during treatment procedures of the breast, lung and thorax. The Butterfly Boards feature a one-piece durable carbon fiber base with two hand posts or an adjustable t-handle.

Bionix Butterfly Boards

- Largest selection of styles to fit all of your immobilization needs
- Available in carbon fiber or acrylic versions
- May be used with SecureVac™ cushions and SecureFoam™
- Attach to any of the Bionix Breast or Torsoboards
- Accommodates Silverman & MoldCare® head cushions
- T-handle grip available on extended versions featuring two axes of adjustability and indexing


To order, or for more information call **Bionix Radiation Therapy** at 1.800.624.6649 or visit us on the web at www.BionixRT.com.



RT-7527
Carbon Fiber
Extended Butterfly Board




RT-7526
Carbon Fiber
Butterfly Board




Butterfly Boards easily attach to
Bionix Breast Boards.

Acrylic Butterfly Boards Include:

RT-5526 Acrylic
Butterfly Board



RT-5527 Extended
Butterfly Board



RT-5528 Narrow Extended
Butterfly Board

