

Max3™ Breast Board

CT Friendly and IMRT ready, the Carbon Fiber **Max3 Narrow Breast Board** is 41cm at the widest point. Minimal loose parts and ergonomic design are highlighted by the fully integrated wedge system and secure Velcro closure straps providing ease of use and transport.



The Max3 Narrow Breast Board features the following benefits:

- Four integrated elevation wedges at 7°, 12°, 17°, and 23°.
- Backstop Cushion to prevent patients from involuntarily sliding downward.
- Able to accommodate patients up to 300lbs.
- Easily accepts Head Only and Head & Shoulder thermoplastic masks.
- Compatible with Silverman head supports.

Breast Board Accessories



Clocking Plate Arm Support
RT-5035-60 Left Support
RT-5035-70 Right Support



RT-5043 Calibrated Arm Cups



*The products you want,
the service you deserve.*

5154 Enterprise Blvd. - Toledo, OH 43612
phone 800.624.6649 www.BionixRT.com



Butterfly Boards



RT-7527 Carbon Fiber Extended Butterfly Board



RT-7526 Carbon Fiber Butterfly Board

All Butterfly Boards include adapter posts to allow for breast board compatibility.
Acrylic Butterfly Boards also available.

SecureFit™ Bars

SecureFit Bars are also available in or order to securely attach the Max3 to your treatment table.



RT-7010 Carbon Fiber SecureFit™ Bar Varian Exact



RT-5010A SecureFit™ Bar Varian Exact



RT-5018 Elekta™ Bar, RT-5016 Siemens Bar



RT-5012A SecureFit™ Bar II Varian Exact

Indexable sliding bar allows for off-center treatments.

*The products you want,
the service you deserve.*

5154 Enterprise Blvd. - Toledo, OH 43612
phone 800.624.6649 www.BionixRT.com

