

# Prone Breast System



## RT-6025 Prone Breast System

The Bionix Prone Breast System is one of the most popular products in the line of comfortable patient positioning devices. The Prone Breast System features a foam cushion set designed to attach to the T-Form Multi-Treatment Base, which can also be used with the Pelvis/Belly PLUS System cushion set.



With each purchase of this system, Bionix is donating a portion of the proceeds to breast cancer research.

Included: 15° Contralateral Breast Wedge Insert



## Accessories



**RT-6125** Box of 50 disposable face covers.



**RT-8000** Multi-Treatment Storage Cart

## BREAST & LUNG POSITIONING



Prone Breast System shown with indexing face cushion.



Features T-Form Multi-Treatment Base; the same base used for the Pelvis/Belly System.

### Features:

- Utilizes prone positioning for breast treatment to deliver less irradiation to the chest wall, lung, and cardiac tissue.
- Available in two height options; 23cm and 15cm.
- Massage-style face cushion supports the head with openings for air flow.
- Reversible, carbon fiber breast bridge holds the contralateral breast comfortably out of the treatment field.
- New hand rest with multiple indexing hand positions.
- Comfortable foam design cradles the patient during treatment cycle.
- Lightweight design; total weight (including T-Form Multi-Treatment Base) is 18 pounds.
- Available as a system with T-Form Multi-Treatment Base, or as a Cushion Set (RT-6025-10), for those customers with a Multi-Treatment Base from the Pelvis/Belly System.
- Locks down onto most treatment couches.

Made in the U.S.A.

**RT-7010** Carbon Fiber SecureFit™ Bar Varian Exact®

**RT-5010A** SecureFit™ Bar Varian Exact®

**RT-5018** Elekta™ Bar, **RT-5016** Siemens Bar

**RT-5012A** SecureFit™ Bar II Varian Exact®  
Indexable sliding bar allows for off-center treatments.