


# SafeStraw™ - A Volume Limiting Drinking Aid



Bionix SafeStraw allows healthcare providers to safely manage swallowed liquid volume and help reduce aspiration of patients with oral motor dysfunction.

## ► PROBLEM:

Patients with oral motor dysfunction have difficulty controlling a liquid bolus and often swallowed prematurely causing them to aspirate.

## ► SOLUTION:

SafeStraw limits volume to a teaspoon, or 6.2 ml's, per sip helping to reduce the likelihood of aspiration by promoting a controllable liquid bolus.

## ► RESULTS:

Limited bolus size helps reduce aspiration risks allowing patients to safely, discreetly, and independently drink from a variety of containers such as cups and bottles.

### PRODUCT OFFERING

**1500 - Box of 12 White, Thin Liquid SafeStraws**

**1510 - Box of 12 Blue, Thick Liquid SafeStraws**

**1520 - Variety Pack** (6 White and 6 Blue SafeStraws)

### POTENTIAL REIMBURSEMENT (U.S. ONLY)

**92526** - For potential reimbursement for treatment of a swallowing dysfunction and/or oral function for feeding.



Scan in the QR code with your web enabled camera phone to view our SafeStraw Video.



©2016 Bionix Development Corp. Patent Pending.



By limiting volume per suck, it has also prevented aspiration in my patients who have decreased oral motor control following a stroke. It makes my patients feel safe and dignified, as they are able to discreetly use the device in restaurants after leaving the hospital.

- Jennifer Hanners, SLP



Obelis s.a • Boulevard Général Wahis 53 • 1030 Brussels, BELGIUM  
Tel: +(32) 2.732.59.54 • Fax: +(32) 2.732.60.03 • E-Mail: mail@obelis.net



Bionix Development Corp. • 5154 Enterprise Blvd. Toledo, OH, 43612, U.S.A.  
Phone: 800-551-7096 • Fax: 800-455-5678 • www.Bionix.com

Phone: 800.551.6810  
www.Bionix.com

**B BIONIX®** GOING-BEYOND